



PDA/FDA Glass Quality Conference

June 4-5, 2012

Renaissance Washington, D.C. Downtown Hotel
Washington, D.C.

**An Alternative Glass Packaging Solution to
Reduce Delamination Risks**

www.pda.org/glass2012

Exhibition: June 4-5 | Courses: June 6-7

Christophe Wagner - SGD



Agenda

- Background:
 - Brief review of Delamination
 - Overview of pharma glass packaging options
- Comparison study between molded glass and current glass packaging solutions
- Observations / Conclusions



Background

SGD

- Independent Molded Glass Producer (formerly Saint-Gobain Desjonquères)
- Dedicated Pharmaceutical glass operations in France and Germany
- R&D lab is located in Mers-Les-Bains Facility, France, where Type I glass is produced.

SGD's evolving view of the delamination issue

- Pre -2010: aware of possible issues with organic acids such as NaHCO₃ 8.4% or Calcium Gluconate
- 2010 ~ mid-2011: Tracking the heightened industry concern about the more severe form of dissolution known as delamination and lamella formation
- Mid-2011 ~ Present: Conducting studies with clients to characterize molded glass chemical durability characteristics; in particular relative to other glass containers they are using

FDA U.S. Food and Drug Administration
Protecting and Promoting *Your* Health

A to Z Index | Follow FDA | FDA Voice Blog

delamination SEARCH

Most Popular Searches

Home Food Drugs Medical Devices Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Radiation-Emitting Products Tobacco Products

Drugs
Home Drugs Drug Safety and Availability

Drug Safety and Availability

- Drug Alerts and Statements
- Importing Prescription Drugs
- Medication Guides

Advisory to Drug Manufacturers: Formation of Glass Lamellae in Certain Injectable Drugs

[3-25-2011] The U.S. Food and Drug Administration (FDA) is announcing the recall of certain injectable drugs containing glass vials. The recall is due to the formation of glass lamellae (glass fragments) in the vials, which have recently been recalled due to the potential for glass fragments to be injected into the patient.

... have been associated with a higher incidence of the formation of glass lamellae:

- ... by tubing process (and thus manufactured under higher heat). These vials are ... glass vials and may shed lamellae more easily.[6] The processing conditions ... glass vials can be designed to mitigate the potential for later delamination.
- ... formulations formulated at high pH (alkaline) and with certain buffers. Common buffers associated with ... information include citrate and tartrate.[7]
- ... the drug product remains exposed to the inner surface of the container. The time duration ... correlation to the potential for glass lamellae formation to occur during the product shelf life.3
- ... with room temperature storage requirements. Drugs stored at room temperature have a ... greater amount of glass lamellae formation than do products stored at colder temperatures.[8]
- ... Termination of the recall has a significant effect on glass stability.4

The referenced literature, below, includes recommended actions to help prevent the formation of glass lamellae. For example, for products "at risk" the vial surface alkalinity can be minimized by proper selection of glass composition (e.g., highly resistant, non-alkaline earth borosilicate glass) appropriate selection and qualification of vendors, and proper quality control of the incoming vials. Accordingly, FDA advises drug manufacturers of products to re-examine their supplier quality management program with the glass vial manufacturers to assure that this phenomenon is not occurring. Further, the Agency reminds finished drug manufacturers to assure that this phenomenon is not occurring. Further, the Agency reminds finished drug manufacturers to assure that this phenomenon is not occurring.

Internet

"... for products "at risk" the vial surface alkalinity can be minimized by proper selection of glass ..."

"Glass vials manufactured ...under higher heat... are less resistant than molded glass vials and may shed lamellae more easily."



Rx-360 is using this flash re

It has been reported... are lit
respect to glass resistance and
are superior to others in pertaining to preventing glass delamination. For example:

“... not all Type I glass is equivalent with respect to glass resistance and delamination... For example: Molded vials are more durable and less susceptible.”

“... increasing surface alkalinity and reducing chemical durability...”

... vial forming makes the glass more durable and less susceptible to delamination.^{3,6} For example:

surface of the vial (where the vial contacts the liquid), thereby increasing surface alkalinity and reducing chemical durability. This physical change continues for the entire life of the vial.⁶

“Listed in order of preference
Best to Worst
1. Molded vial ...”



pH		>7.0	5
Buffer Type	N/A	Citrate	4, 8
Ionic Strength	N/A	>100 mM NaCl	8
Listed in order of preference Best to Worst			
Vial Configuration	1. Molded vial 2. Silica coated tubing vial 3. Regular tubing vial 4. Ammonium sulfate treated tubing vial		4,5,6,7,8




Overview of pharma glass packaging options

Tubing Glass

- 2 step process:
 - Cane manufacturing
 - Converting
- Capabilities:
 - Vials
 - Cartridges
 - Syringes



Molded Glass 1 step process



SGD Capabilities:

- Vials and IV bottles from 3 ml to 1 L. Neck finish 20 mm and higher
- Can produce non round vials and bottles

Comparison study between molded glass and current glass packaging solutions For Type I Vials 5 ml and 10 ml

Study overview

Objective:	To characterize leachable/extractable profile of molded and tubing vials across a range of pH; and to further our understanding of the referenced correlations to hydrolytic resistance		
Scope:	Sample id	Description	
	M-5	Asolvex® molded (SGD) vials;	5ml x 20mm
	M-10	Asolvex® Molded (SGD) vials;	10ml x 20mm
	T-5	Tubing vials, European Supplier 1;	5ml x 20mm
	T-10	Tubing vials, European Supplier 2;	10ml x 20mm

Study overview (cont.)

Analysis:	Methods
1- Mass Composition	X-Ray Fluorescence Spectrometry
2- Surface Composition	Time-of-Flight Secondary Ion Mass Spectrometry (TOF-SIMS)
3- Hydrolytic Resistance	Grain and Surface per European Pharmacopeia, 3.2.1
4- Extractables	Aqueous extraction at elevated temperature and across a pH range of 1 – 10 + Emission Spectrometry ICP

Study 1: Mass Composition Analysis

- Method: X-Ray Fluorescence Spectrometry
 - Vials are cut in pieces
 - Samples flattened at 750°C
 - Surface is polished
 - X-Ray Fluorescence on 34mm diameter samples
 - FX S8 TIGER BRUKER

Reminder: Type I glass composition

- NEUTRAL GLASS is an alkaline borosilicate glass with main components of (typical moulded glass composition):
 - Network Formers :

$\text{SiO}_2 + \text{Al}_2\text{O}_3$	-	73%
B_2O_3	-	12%
 - Network Modifiers:

$\text{Na}_2\text{O}; \text{K}_2\text{O}$	-	10%
$\text{CaO}; \text{BaO}; \text{ZnO}$	-	5%
- NEUTRAL GLASS may be composed of 2 primary phases
 1. Silica-rich phase with low alkaline content
 2. Boron-rich phase with most alkaline elements of the glass; it may be separated into micro-droplets within the silica rich matrix, depending on the composition

Results / Conclusion

(%)	Molded	Tubing 1	Tubing 2
Network Formers	85.7	90.2	91.1
Network Modifiers	14.2	9.6	8.7

- Stronger network for bulk tubing glass, less modifiers
- Network modifiers needed to soften the glass to shape the vials for molded glass

Main elements (%)	Moulded Flint	5ml Tubing 1	10ml Tubing 2
SiO ₂	69,1	70,8	74,3
Na ₂ O	6,1	7,1	7,2
K ₂ O	3,1	1,2	0,0
CaO	1,1	1,2	1,5
MgO	0,0	0,2	0,0
Al ₂ O ₃	4,0	7,3	5,6
Fe ₂ O ₃	0,02	0,03	0,02
B ₂ O ₃	12,6	12,1	11,2
BaO	2,8	0,1	0,0
TiO ₂	0,02	0,01	0,03
ZnO	1,1	0,0	0,0

Study 2: Surface Composition Analysis - SIMS

- Surface SIMS analysis by Time-of-Flight Secondary Ion Mass Spectrometry (ToF-SIMS)
 - 4 glass vial samples : 2 molded and 2 tubing vials
 - ToF SIMS Profile by alternating analysis and abrasion cycles
 - Analysis:
 - Primary Ions Bi_1^+ 25 keV, $I = 1\text{pA}$
 - Surface analyzed $100 \times 100 \mu\text{m}^2$, 128x128pixel
 - Positive Secondary Ions analyzed
 - Abrasion:
 - Primary Ions O_2^+ 500eV, $I = 100\text{nA}$
 - Surface : $300 \times 300 \mu\text{m}^2$
 - Cycle
 - Analysis : acquisition of 1 scan (time of max flight = $100 \mu\text{s}$)
 - Abrasion : 1.6s, Pause : 1s

Glass Composition : from surface to internal




B+

Na+

TUBING 5ml

Al+

K+







B+

Na+


MOLDED 5ml

Al+

K+



Glass Composition : from surface to internal



Conclusion – Surface Composition

- All samples show a different surface composition
- Small and curved samples may explain different bulk compositions between the bottom and the side wall
- More surface composition differences between side wall and bottom for tubing vials
- Sodium depletion at surface on the vial bottom for tubing
- Sodium depletion during forming for Asolvex Type I glass, both on bottom and on side walls (blowing effect)

Study 3 : Hydrolytic Resistance Testing

- Standard test for Pharma Glass - Hydrolytic stability, expressed by the resistance to the release of soluble mineral substances into water under the prescribed conditions of contact between :
 - the inner surface of the container (Test A, surface test according to European Pharmacopeia, 3.2.1)
 - glass grains and water (Test B, glass grain test according to European Pharmacopeia, 3.2.1)
- The hydrolytic resistance is evaluated by titrating released alkali.
- The glass grain test is performed on crushed glass pieces, so represents the chemical resistance of the bulk glass

Hydrolytic Resistance Comparison in (ml) HCl N/100

	Type I Molded	Tubing T-5
Grain Hydrolytic Resistance (ml)	0.53	0.43

- Better grain resistance for Tubing than molded because more network formers and less modifiers, Type I Limit 1 ml


	Type I Molded M-5	Tubing T-5	Type I Molded M-10	Tubing T-10
Vol 90% (ml)	8.1	8.3	12.25	12.4
Type I Limit	1	1	0.8	0.8
Surface Hydrolytic Resistance (ml)	0.15	0.50	0.17	0.41

- More critical for product interaction
- All vials are under type I surface limit
- Better surface Hydrolytic resistance for molded vials

Study 4: Extractables evaluation


- Autoclave solution analysis with ICP
- **Solution Preparation**
 - Deionized water pH (18 MΩ.cm resistivity) adjusted :
 - with HCl for acid pH
 - with NaOH for base pH
- **Vials Extraction**
 - filled at nominal capacity with the solution
 - Vials in autoclave at 121°C for 1h, Eur. Pharma. HR cycle , 3 to 5 for each pH
- **ICP Preparation**
 - Acidification HNO₃ Suprapur 2% before ICP measurement
 - Equipment Calibration with certified PE multielements solution and acidification HNO₃ Suprapur 2%
- **Results**
 - Equipment : Emission Spectrometry ICP (Perkin Elmer Optima 7300 DV)
 - The blank solution is analyzed and subtracted from the autoclaved solutions.

Vial comparison : Total Extractables – 5ml



- Less elements extracted with Molded vials, for all pH
- Higher pH (10 or more) causes higher extractions

Vial comparison : Total Extractables – 10ml



- Higher pH (10 or more) causes higher extractions
- Less extraction in volume for bigger vial, surface/volume ratio lower

Extractables Analysis by element– 5ml

	PH=1		PH=2		PH=4		PH=6		PH=7		PH=8		PH=9		PH=10	
Extracted Elements (µg/L)	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5	Moulded M-5	Tubing T-5
Si	241	1632	203	1320	188	1118	368	3216	640	3443	818	3253	1079	3447	3481	6315
Na	272	1913	246	1647	185	1162	137	883	146	800	158	881	209	1026	471	1735
K	126	213	111	190	75	119	65	111	78	99	87	127	109	135	334	250
Ca	136	326	74	281	60	122	94	218	93	199	131	263	143	381	229	606
Mg	6	7	2	7	2	4	3	6	4	6	3	9	5	10	5	15
Al	58	771	54	541	34	221	62	509	4	84	87	511	140	630	339	1068
Fe	20	5	-	4	-	2	8	3	0	0	7	3	13	10	10	14
B	123	1058	99	939	51	691	62	639	73	585	99	578	158	675	421	1075
Ba	64	4	52	22	34	15	39	21	32	9	68	47	89	42	228	34
Ti	1	1	1	1	0	0	1	0	0	0	0	0	1	2	2	2
Zn	58	3	45	15	34	7	37	14	34	4	35	15	69	21	130	25
Extractables Total (µg/L)	1 105	5 931	887	4 967	663	3 459	872	5 618	1 101	5 227	1 491	5 684	2 012	6 377	5 648	11 137

- No visible attack of the glass, no flake (methylene blue test shows nothing)
- Different local / surface glass compositions with tubing may cause higher extractions

Extractables Analysis by element– 10ml

	PH=1		PH=2		PH=4		PH=6		PH=7		PH=8		PH=9		PH=10	
Extracted Elements (µg/L)	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10	Moulded M-10	Tubing T-10
Si	201	781	141	493	146	275	310	1535	312	1575	465	2116	497	1920	2320	5688
Na	200	1492	199	1195	213	752	93	567	117	540	112	741	129	754	329	1781
K	96	37	101	8	119	5	41	0	68	27	63	8	81	7	208	10
Ca	39	272	47	165	66	56	37	141	28	109	38	194	46	219	81	365
Mg	1	5	2	2	3	0	1	2	1	2	2	2	2	4	2	5
Al	49	457	49	305	33	68	46	277	1	68	61	368	81	396	244	964
Fe	6	11	6	7	4	2	6	11	0	0	3	5	3	6	3	13
B	86	707	83	506	33	212	51	341	41	304	63	447	69	431	252	978
Ba	44	0	53	0	29	0	30	4	14	0	44	0	45	6	141	0
Ti	0	2	0	1	0	0	0	0	0	0	0	1	1	1	1	5
Zn	24	15	29	14	18	7	19	15	14	8	22	13	28	21	67	21
Extractables Total (µg/L)	742	3 777	709	2 694	662	1 374	631	2 890	592	2 629	870	3 894	979	3 763	3 644	9 829

- Values lower than 20 µg/L may not be significant (pollution ?),
- Still refining our test method, ICP detection limit on the blank solution $3\sigma < 4\mu\text{g/L}$ (σ calculated on 10 measurements of the blank solution)
- Vial to vial variation +/- 10%

Comments on Extractables

- **Tubing** : more Na and Ca extracted, but also Al, Si and B which are the glass network formers
- **Molded** : more K (not in the tubing 10ml glass composition) and Ba (traces in the tubing glass composition), which are mainly glass modifiers and less impacting the glass chemical robustness
- Bulk hydrolytic resistance is good for tubing, but surface resistance is not at the same level
- Local changes in glass compositions (processing effect) may explain some of the increased extraction

Conclusions

- Delamination is a last stage indicator of heavy extraction of glass by the product
- Product interaction with the glass depends on glass composition AND how it was formed
- Results seem to indicate 1 step forming of molded seems to extract glass formers less readily than 2 step tubing process
- Tubing glass starts off better at cane stage but chemical robustness is impacted by converting step, which can differ from 1 supplier to another
- Due to its chemical robustness, molded could be considered as an alternative in high aggressive extraction conditions
- Further comparison studies with real products are being conducted with clients to confirm this theory

Thank you for your attention !

Acknowledgements

Work done by SGD lab in Mers-les-Bains, France

- Didier Pichard
- Joel Bourjot
- Sébastien Dussardier

SIMS profiles done by CRITT